# Liqui-Phil<sup>TM</sup> General Manual

Copyright © 2016

#### WARNING

Do not allow filled phantoms to be exposed to heat or freezing. This product is not compatible with some chemicals. Review chemicals to be used in the phantom for compatibility with CAB and polycarbonate materials prior to use. It is advised not to expose this product to open flame or high temperature (over 38° Celsius or 100° Fahrenheit) heating elements.

#### WARRANTY

THE PHANTOM LABORATORY INCORPORATED ("Seller") warrants that this product shall remain in good working order and free of all material defects for a period of one (1) year following the date of purchase. If, prior to the expiration of the one (1) year warranty period, the product becomes defective, Buyer shall return the product to the Seller at:

By Truck The Phantom Laboratory Incorporated 2727 State Route 29 Greenwich, NY 12834

Or By Mail PO Box 511 Salem, NY 12865-0511

Seller shall, at Seller's sole option, repair or replace the defective product. The Warranty does not cover damage to the product resulting from accident or misuse.

IF THE PRODUCT IS NOT IN GOOD WORKING ORDER AS WARRANTED, THE SOLE AND EXCLUSIVE REMEDY SHALL BE REPAIR OR REPLACEMENT, AT SELLER'S OPTION. IN NO EVENT SHALL SELLER BE LIABLE FOR ANY DAMAGES IN EXCESS OF THE PURCHASE PRICE OF THE PRODUCT. THIS LIMITATION APPLIES TO DAMAGES OF ANY KIND, INCLUDING, BUT NOT LIMITED TO, DIRECT OR INDIRECT DAMAGES, LOST PROFITS, OR OTHER SPECIAL, INCIDENTAL, OR CONSEQUENTIAL DAMAGES, WHETHER FOR BREACH OF CONTRACT, TORT OR OTHERWISE, OR WHETHER ARISING OUT OF THE USE OF OR INABILITY TO USE THE PRODUCT. ALL OTHER EXPRESS OR IMPLIED WARRANTIES, INCLUDING, BUT NOT LIMITED TO, THE IMPLIED WARRANTY OF MERCHANT ABILITY AND FITNESS FOR PARTICULAR PURPOSE, ARE HEREBY DISCLAIMED.

## The Phantom Laboratory

### Contents

Warning 1

Introduction 4

Care and Maintenance 5

Chemical Compatibility Chart 6

Optional Accessories 12

-Tumor Rotation Ball 12

-Optional Tumor Vessels 12

#### Introduction

The Phantom Laboratory designs and manufactures its own phantoms as well as custom and OEM phantoms. These phantoms are used for numerous research activities where scientists need to hold liquid solutions. Because each application is different, and each medical facility has its own unique set of requirements, The Phantom Laboratory does not set forth operational protocols, nor do we make specific recommendations on the frequency and expected results of evaluations. Review the local governing regulations, the needs of your radiologists and physicists, and your system manufacturer's recommendations when developing your measurement protocols. In order to prevent damage to your phantom, review the chemical compatibility chart that begins on page 6 in order to determine which chemicals can be used in these phantoms.

If you have questions about the phantom, please contact The Phantom Laboratory at the following address:

The Phantom Laboratory, Incorporated PO Box 511, Salem NY 12865 USA Phone: 800-525-1190 or 518-692-1190

Fax: 518-692-3329

Email:info@phantomlab.com

#### **Care and Maintenance**

The Liqui-Phil™ Phantom shell is fabricated from 0.125-inch-thick cellulous acetate butyrate (CAB), a transparent plastic chosen for its strength and low water absorption. The shell is mounted on durable polycarbonate end plates. Both of these materials are susceptible to solvents, therefore we recommend the use of mild soap and water when cleaning. Please note that water temperatures above 100° Fahrenheit should not be used as it will cause the plastic material to become deformed, thus destroying the phantom.

The Liqui-Phil™ Phantom is designed to be filled with water. Please refer to the Chemical Compatibility Chart on pages 6-11 of this manual prior to use. Some chemicals will destroy the Phantom's shell and, if used, will render the warranty invalid.

You should never allow a filled phantom to freeze. If water is left in the phantom for an extended period of time, bacteria and other microorganisms may begin to grow and cause the water to become cloudy. If clouding occurs, the phantom should be emptied and washed with a mild detergent and water. Before storing the phantom, empty all water from it. To prevent mold build-up, the phantom should be stored empty, clean, and dry with the port plugs removed.

Cellulose acetate butyrate will yellow when exposed to ultraviolet light. This discoloration will not affect the phantom's performance.

The Liqui-Phil™ Phantom may contain o-rings. To ensure easy function and a tight seal, always lubricate o-rings with petroleum jelly or silicon lubricant periodically during use and after each cleaning.

## **Chemical Compatibility Chart**

	Perc	ent Inc	rease	Observed Condition
Reagent	Time Exposed Weight			of Plastic
Chemicals	<b>F</b>			
Acids				
Acetic, 5%	1 year	3.6	2.1	Slightly softened
Acetic, 10%	2 months	5.2	2.4	Slightly softened
•Acetic, 30%	2 months	13.6	8.6	Softened and swollen
Boric, 5%	2 days	1.3	0.0	Unchanged
•	8 days (38°C, 100°F)	2.0	0.1	Slightly stained
	4 months (60°C, 140°F)	1.6	0.8	Slightly softened
Citric, 10%	1 year	1.4	0.6	Unchanged
	4 months (60°C, 140°F)	_	_	Surface attacked
Fluosilicic, 10%		4.5	1.2	Unchanged
Fluosilicic, 28%		4.7	3.6	Unchanged
Formic, 3%	20 days	_	_	Unchanged
Hydrochloric, 1		0.9	0.5	Surface slightly attacked
•Hydrofluoric,		10.3	5.5	Slightly swollen and softened
•Hydrofluoric,		_	_	Dissolved
Lactic, 50%	2 days	1.6	0.5	Unchanged
•Nitric, 10%	8 months	_	_	Decomposed
Oleic	1 year	2.3	1.5	Unchanged
Phosphoric, 309	•	1.3	0.8	Unchanged
Phosphoric, 509		1.6	0.8	Unchanged
•Phosphoric, 75		_	_	Partially decomposed
Pyrogallic, 4%	1 week	2.6	1.1	Stained yellow
Stearic	1 week	_	_	Unchanged
Sulfuric, 3%	1 year	1.6	1.0	Slightly discolored
Sulfuric, 10%	1 year	1.5	0.7	Slightly discolored
•Sulfuric, 20%	1 year	0.9	0.3	Slightly softened, surface
attacked •Sulfu	•	0.4	0.3	Surface attacked
•Sulfuric, 94%	.,	_	_	Disintegrated
	4 months (38°C, 100°F)	2.8	1.2	Unchanged
Trichloroacetic		3.3	0.5	Unchanged
•Trichloroaceti		9.3	3.1	Softened
Alcohols, Mono	· ·			
n-Amyl	2 days	3.1	3.0	Unchanged
•tert-Amyl	2 days	14.0	11.3	Softened, tacky
•n-Butyl	2 days	6.5	7.2	Swollen
•sec-Butyl	2 days	7.2	10.7	Swollen
•tert-Butyl	2 days	3.6	3.3	Slightly softened
•Diacetone	•	_	_	Dissolved
•Ethyl (denatur	red) 2 days	23.0	24.7	Softened
•Ethyl, 50%	1 week	13.4	11.6	Softened
•2-Ethylhexyl	1 week	_	_	Swollen
Hydrocarbons				
Gas, Natural, ar	romatic-free 1 year	_	_	Showed slight decrease in
				tensile strength and
				increase in impact strength
Gas, Natural, 59	% aromatic content 23	days	_	Showed slight decrease in
				tensile strength and
				increase in impact strength
				·

<sup>•</sup>Indicates that material is generally unsatisfactory for use in contact with Tenite butyrate under the conditions of this test

Chemical

**Compatibility Chart continued** 

	Percen	t Incre	ase	Observed Condition	
Reagent	Time Exposed Weight			of Plastic	
Chemicals					
•Isoamyl	$2  ext{ days}$	2.0	2.1	Very slightly softened	
•Isopropyl	2 days	23.4	<b>25.1</b>	Softened, tacky	
•Methyl		_	_	Dissolved	
Methyl, 5%	1 year	2.0	1.2	Slightly softened	
•n-Propyl	2 days	<b>15.0</b>	4.4	Slightly softened	
•Tetrahydrofu	rturyi dric and Trihydric	_	_	Dissolved	
•Diethylene G		8.2	6.1	Softened	
	ediol-1,3 2 days (38°C, 1			Unchanged	
Ethylene Glyc		4.2	2.1	Unchanged	
Glycerin	1 year	0.0	0.4	Unchanged	
Propylene Gly		0.4	0.0	Unchanged	
•Triethylene		8.6	6.7	Softened	
Bases					
•Ammonium I	Hydroxide, 10%				
	2 months	21.9	12.9	Softened	
Calcium Hydr	oxide, saturated solutio				
~	1 week	0.7	0.7	Unchanged	
-	oxide, 1% 1 year	1.0	0.6	Unchanged	
	roxide, 10% 8 months	3.2	2.2	Brittle	
Trimethylben	zyl Ammonium Hydroxi	,	0.0	Unchanged	
Esters	17 days	1.1	0.0	Unchanged	
•n-Butyl Aceta	ate	_	_	Dissolved	
•sec-Butyl Ace		_	_	Dissolved	
Di-2-Ethylhex	0.9	0.3	Unchanged		
Di-2-Ethylhex	0°C, 122	P°F)	Small gain		
•Ethyl Acetate			_	Dissolved	
•Ethyl Lactate		_	_	Dissolved	
•Ethyl Propio	_	_	Dissolved		
•Ethylene Gly		_	Dissolved		
	col Monomethyl Ether	Acetate	_	Dissolved	
•Isoamyl Acet		_	_	Dissolved	
•Isobutyl Acet		_	_	Dissolved	
•Isopropyl Ac		_	_	Dissolved	
•Methyl Aceta		_	_	Dissolved	
•n-Propyl Ace	tate	_	_	Dissolved	
Ethers • Dichloro Diethyl Ether — Dissolved					
•Diethyl Ethe		46.0	50.0	Considerably swollen	
Di-Isopropyl I		0.8	1.1	Unchanged	
Ether-Alcohol	•			<b>g</b>	
•Ethylene Gly	col Monoethyl Ether	_	_	Dissolved	
	col Monomethyl Ether	_	_	Dissolved	
<ul><li>Toluene</li></ul>	2 days	39.3	<b>54.9</b>	Softened	
•Xylene	1 week	41.5	33.2	Softened	
Hydrocarbons, Halogenated					
•Carbon Tetrachloride 2 days		14.8	6.8	Surface slightly softened	
•Chlorobenzene		_	_	Dissolved	
• Chlorobromomethane		_	_	Dissolved	
• Chloroform		_	_	Dissolved	
	enzene 3 days	_	_	Softened and swollen	
		 11.1		Swollen	
1				SWOHOH	
Chemical Compatibility Chart continued					

		Perce	ent Incre	ease	Observed Condition
Reagent	Time Ex	posed Weigl	ht Thick	ness	of Plastic
Chemicals					
•Ethylene Chlo			_	_	Dissolved
•Methylene Ch			_	_	Dissolved
•Propylene Ch		_	_	_	Dissolved
•s-Tetrabromo		3 days	_	_	Softened, swollen, and tacky
•Tetrachloroet			_	_	Dissolved
•Tetrachloroet		12 days	_	_	Badly swollen
•Trichloroethy	lene	1 day	_	_	Badly swollen
Ketones					D' 1 1
•Acetone			_	_	Dissolved
•Cyclohexanor			_	_	Dissolved
•Di-Isopropyl			_	_	Dissolved
•Methyl Ethyl			_	_	Dissolved
•Methyl n-Buty			_	_	Dissolved
•Methyl Isobut	tyl Keton	ie	_	_	Dissolved
•Phorone			_	_	Dissolved
Salts	4-4- <b>D</b>		1		
Aluminum Ace		ac, 33% watei 2 months	_	0.0	II. ah au mad
Al Ch1	-		1.8	0.6	Unchanged
Aluminum Chl			1.5	0.7	Unchanged Aluminum
Chloride, satur		ution 2 months	Λ1	0.0	II. ah an mad
A1 C14			0.1	0.0	Unchanged
Aluminum Sulf			1.7	1.1	Unchanged
Ammonium Bit		saturated so 1 month	2.3		Slightly blooghed
Ammonium Ch				_	Slightly bleached
Ammonium Ch		aturated som 1 month	2.1	0.8	Unchanged Ammonium
Nitrate, solid		1 month 1 week	0.2	0.8	Unchanged Ammonium
Nitrate, sond Nitrate, 10%		1 week 1 week	1.7	0.2	Unchanged Ammonium
Sulfate, solid		1 week 1 year	0.1	0.6	Unchanged
Ammonium Su		•	1.3	0.5	Unchanged
Calcium Chlor			1.5	0.9	Unchanged
Calcium Chlor		-	0.4	0.0	Unchanged
Calcium Hypod			6.0	-3.8	Softened and cracked
Calcium Hypod			0.8	0.0	Unchanged
Calcium Phosp				0.0	Ononungou
culcium i nosp		1 year	1.7	0.5	Unchanged
Calcium Phosp		•		0.0	0.1011411904
ошения и поор		1 year	-0.6	0.6	Unchanged
Calcium Phosphate, Tribasic, solid					
		1 year	-0.6	0.6	Unchanged
Heptane		1 year	1.6	2.5	Unchanged
Hexane		1 week	_	_	Unchanged
Propane, gas	9	2 months	0.3	0.6	Unchanged
Propane, liquio	d :	2 months	1.4	4.6	Unchanged
Calcium Sulfat		ım), solid			8
		1 year	-0.1	0.6	Unchanged
Copper Sulfate			ths 1.7	0.6	Unchanged
Copper Sulfate	, saturat	ed solution			
		2 months	1.7	0.9	Unchanged
<b>Cuprous Chlor</b>	ide, solid	l 1 week	1.5	_	Unchanged
Ferric Ammoni					
1 week (38°C, 100°F, 0.3 — Unchanged					
80% RH) Ferric			hs 2.0	0.8	Unchanged
Ferric Chloride		2 months	1.7	1.0	Unchanged
Ferric Chloride	e, <b>40</b> %	$2~\mathrm{months}$	1.3	0.4	Unchanged
<b>Chemical Con</b>	npatibil	ity Chart co	ontinue	d	

	nt Increa		<b>Observed Condition</b>		
Reagent Time Exposed Weight			of Plastic		
Ferric Chloride, saturated solution 1	month (		Unchanged		
Lithium Bromide, solid 1 week	- 0.7	0.0	Unchanged		
Lithium Bromide, 50% 1 week	0.0	0.0	Unchanged Magnesium		
Carbonate, 2.5% 2 days	1.6	1.0	Unchanged		
Potassium Aluminum Sulfate (alum),		0.8	Unchanged		
4 months (38°C, 100°F) Potassium Bromide, 3%	1.9	0.0	Unchanged		
3 days (38°C, 100°F)	1.3	_	Unchanged		
Potassium Chloride, solid 1 year	0.1	0.5	Unchanged		
Potassium Chloride, 10% 1 year	1.7	0.4	Unchanged Potassium Chrome		
Alum, 10% 3 days (38°C, 100°F) 1.3	_	Unch	anged		
Potassium Cyanide, 10% 2 months	1.4	0.3	Slightly discolored (brown)		
Potassium Cyanide, saturated solution	n				
2 months	0.5	0.0	Slightly discolored (brown)		
Potassium Ferricyanide, 10% 4 days	_	_	Unchanged		
Potassium Sulfate, solid 1 year	0.1	0.6	Unchanged		
Potassium Sulfate, 10% 1 year	1.4	0.4	Unchanged		
Silver Nitrate, 2.5% 2 days	1.5	0.0	Unchanged		
Sodium Acetate, 3% 3 days (38°C, 100°			Unchanged Unchanged		
Sodium Aluminum Sulfate, solid 1 w Sodium Bicarbonate, 2.5% 2 days	1.7	$\begin{array}{c} 0.4 \\ 0.5 \end{array}$	Unchanged		
Sodium Bisulfate, solid	1.7	0.0	Chenangeu		
1 week (38°C, 100°F,	0.1	_	Unchanged		
80% RH) Sodium Bisulfate, 1%	0.1		Chonangea		
3 days (38°C, 100°F)	1.3	_	Unchanged		
Sodium Bisulfite, 20% 1 week	2.1	0.8	Unchanged		
Sodium Borate, 2.5% 2 days	1.5	0.5	Unchanged		
Sodium Carbonate, solid			_		
1 week (38°C, 100°F,	3.9	_	Unchanged		
80% RH) Sodium Carbonate, 2.5% 1 year 1.3 0.9			Unchanged Sodium		
Carbonate, 6% 3 days (38°C, 100°F	1.2	_	Unchanged		
Sodium Carbonate, 10% 1 year	_	_	Unchanged		
Sodium Chloride, 2.5% 1 year	_	_	Unchanged		
Sodium Chloride, 10% 1 year	1.3	0.5	Unchanged		
Sodium Chloride, saturated solution					
2 months	0.8	0.3	Unchanged		
Sodium Chloride, saturated solut	ion				
2 months (60°C, 140°F) 0.9 0.9 Unchanged					
Sodium Chromate, saturated solution					
1 week	0.6	0.2	Unchanged		
Sodium Cyanide, 10% 2 months	1.0	0.3	Unchanged		
Sodium Cyanide, saturated solution					
		- 0.4	Unchanged		
Sodium Ferrocyanide, solid 1 wee		_	Unchanged		
Sodium Fluoride, 4% 1 month	2.5	_	Unchanged		
Sodium Hypochlorite, 30% 13 day		- 2.1	Unchanged		
Sodium Nitrate, solid 2 months	0.1	0.1	Unchanged		
Sodium Nitrate, 10% 1 year	1.2	0.4	Unchanged		
Sodium Nitrate, saturated solution					
2 months	0.9	0.4	Unchanged		
		0.4	Unchanged Unchanged		
,					
Sodium Silicate, saturated solution					
2 months	1.2	0.1	Unchanged		
Chemical Compatibility Chart continued					

Reagent   Time Exposed   Weight   Thickness   Sodium Sulfite, 10%   1 week   2.1   0.9   Unchanged   Sodium Thiosulfate, 24%   3 days (38°C, 100°F)   1.2   —   Unchanged   Time thy Boundary   1 month (50°C, 122°F)-0.7   0.1   Unchanged   Unchanged   Time thy Benzyl Ammonium Chloride, 5%   17 days   1.1   0.1   Unchanged   Unch		Percen	t Increa	ase (	Observed Condition	
Sodium Thiosulfate, 20% 13 days   1.1   0.0   Unchanged	Reagent Time	Exposed Weight	Thickr			
Sodium Thiosulfate, 20% 13 days   1.1   0.0   Unchanged	_	_			Unchanged	
Sodium Thiosulfate, 24%   3 days (38°C, 100°F)   1.2   — Unchanged			1.1	0.0		
Tetra (2-Ethylbutyl) Silicate					S .	
Tetra (2-Ethylbutyl) Silicate		·	1.2	_	Unchanged	
Trimethyl Benzyl Ammonium Chloride, 5%	-	· ·			<u> </u>	
17 days			F)-0.7	- 0.1	Unchanged	
Zinc Chloride (hydrous salt) 1 week 0.5	Trimethyl Benzyl A	mmonium Chl	oride, 8	5%	_	
Zinc Chloride, saturated solution		17 days	1.1	0.1	Unchanged	
1 week	Zinc Chloride (hyd	rous salt) 1 we	e <b>k 0.5</b>	0.0	Unchanged	
Zinc Oxide, solid	Zinc Chloride, satu	rated solution				
Miscellaneous Chemicals and Gases Ammoniated Mercury 1 week (60°C, 140°F)—  *Aniline — — — — — — — — — — — — — — — — — — —		1 week	1.4	0.8	Slightly etched	
Ammoniated Mercury 1 week (60°C, 140°F)—  *Aniline — — — Dissolved  *Benzaldehyde — — — Dissolved  *Butadiene-1,3, liquid 6 months 19.3 26.4  *Butadiene-1,3, gas 1 month 2.7 2.3 Unchanged  *Carbon Disulfide 1 week 25.8 1.6 Softened and swollen  *Carbon Disulfide, saturated atmosphere  2 days 17.4 11.8 Warped  *Chlorine, dry 1 week 8.8 2.2 Crazed and brittle  *Chlorine, moist 1 week 7.8 0.1 Crazed and brittle  *Chlorine, saturated solution 1 week — — Dissolved  *I.,4-Dioxane — — Dissolved  *Ethylene Oxide, gas  10 minutes (41°C, 105°F)——  *Ethylene Oxide, gas 1 day 20.9 25.6 Swollen and softened  Formaldehyde, 4% 10 min per day 0.2 — Unchanged  *Formaldehyde, 35% 2 months 13.0 6.7  *Furfural — — Dissolved  Hydrogen Peroxide, 5% 2 days 1.4 1.3 Unchanged  Hydrogen Peroxide, 5% 2 days 1.4 1.3 Unchanged  Hydrogen Sulfide, moist 2 months 2.3 0.8 Hydrogen Sulfide, saturated solution  2 months 5.9 2.3 Unchanged  Hydroquinone, 20 g per gallon 1 week 2.4 1.0 Sightly stained yellow  *Methyl Methacrylate Monomer — — Dissolved  Unchanged  Hydroquinone, 0.05-0.15 ppm 45 days (outdoors)— — Dissolved  Ozone, 0.05-0.15 ppm 45 days (outdoors)— — Dissolved  *Styrene Monomer — — Dissolved  *Sulfur Dioxide, dry 2 months 19.4 8.6  *Sulfur Dioxide, moist 2 months 31.9 10.2  *Sulfur Dioxide, saturated 2 months 23.2 18.1  Unchanged  *Sulfur Dioxide, saturated 2 months 23.2 18.1	Zinc Oxide, solid	1 week	_	_	Unchanged	
*Aniline  *Benzaldehyde  *Butadiene-1,3, liquid 6 months 19.3 26.4  Butadiene-1,3, gas 1 month 2.7 2.3  *Carbon Disulfide 1 week 25.8 1.6  *Carbon Disulfide, saturated atmosphere  2 days 17.4 11.8  *Chlorine, dry 1 week 8.8 2.2  *Chlorine, moist 1 week 7.8 0.1  *Chlorine, saturated solution 1 week — Considerably softened and swollen  *In minutes (41°C, 105°F) — Unchanged  *In minutes (41°C, 105°F) — Unchanged  *Swollen and softened brittle  *Chlorine, saturated solution 1 week — Dissolved  *In minutes (41°C, 105°F) — Unchanged  *Formaldehyde, 35% 2 months 13.0 6.7  *Formaldehyde, 35% 2 months 13.0 6.7  *Furfural — Dissolved  Hydrogen Peroxide, 5% 2 days 1.4 1.3  Hydrogen Sulfide, dry 2 months 2.3 0.8  Hydrogen Sulfide, moist 2 months 3.0 1.7  Hydrogen Sulfide, saturated solution 2 months 5.9 2.3  Hydroquinone, 20 g per gallon 1 week 2.4 1.0  *Methyl Methacrylate Monomer — Dissolved  Hydroquinone, 20 g per gallon 1 week 2.4 1.0  *Methyl Methacrylate Monomer — Dissolved  Dissolved  Dissolved  Hydroquinone, 20 g per gallon 1 week 2.4 1.0  *Methyl Methacrylate Monomer — Dissolved  Ozone, 0.05-0.15 ppm 45 days (outdoors)— Unchanged  *Sulfur Dioxide, dry 2 months 19.4 8.6  *Sulfur Dioxide, dry 2 months 19.4 8.6  *Sulfur Dioxide, saturated 2 months 23.2 18.1  Dissolved  Swollen and softened  Unchanged  Swollen and softened  Unchanged  Vunchanged  Vunchanged  Unchanged  Unchanged  Unchanged  Unchanged  Unchanged  Unchanged  Vellowed  Dissolved  Unchanged  Vellowed  Phenol 1 week — Decomposed  Styrene Monomer — Dissolved  Unchanged  Swollen, slightly warped  *Sulfur Dioxide, moist 2 months 31.9 10.2  *Sulfur Dioxide, saturated 2 months 23.2 18.1	Miscellaneous Cher	micals and Gas	es			
*Benzaldehyde *Butadiene-1,3, liquid 6 months 19.3 26.4 Swollen and softened Butadiene-1,3, gas 1 month 2.7 2.3 Unchanged *Carbon Disulfide 1 week 25.8 1.6 Softened and swollen  *Carbon Disulfide, saturated atmosphere  2 days 17.4 11.8 Warped  *Chlorine, dry 1 week 8.8 2.2 Crazed and brittle  *Chlorine, moist 1 week 7.8 0.1 Crazed and brittle  *Chlorine, saturated solution 1 week — Considerably softened and swollen  **I,4-Dioxane — Dissolved  *Ethylene Oxide, gas  10 minutes (41°C, 105°F) — Unchanged  *Ethylene Oxide, gas 1 day 20.9 25.6 Swollen and softened  *Formaldehyde, 4% 10 min per day 0.2 — Worlden 13.0 6.7 Swollen and softened  *Formaldehyde, 35% 2 months 13.0 6.7 Swollen and softened  *Hydrogen Peroxide, 5% 2 days 1.4 1.3 Unchanged Hydrogen Sulfide, dry 2 months 2.3 0.8 Unchanged  Hydrogen Sulfide, moist 2 months 3.0 1.7 Unchanged  Hydrogen Sulfide, saturated solution  2 months 5.9 2.3 Unchanged  Hydroquinone, 20 g per gallon 1 week 2.4 1.0 Slightly stained yellow  *Methyl Methacrylate Monomer — Dissolved  Dissolved  Ozone, 0.05-0.15 ppm 45 days (outdoors)— Unchanged  *Sulfur Dioxide, dry 2 months 19.4 8.6 Swollen, slightly warped  *Sulfur Dioxide, moist 2 months 31.9 10.2 Considerably swollen and warped  *Sulfur Dioxide, saturated 2 months 23.2 18.1 Swollen and warped	Ammoniated Mercu	ıry 1 week (60°	C, 140°	<b>F</b> )——		
*Butadiene-1,3, liquid 6 months			_	_		
Butadiene-1,3, gas 1 month 2.7 2.3 Unchanged  *Carbon Disulfide 1 week 25.8 1.6 Softened and swollen  *Carbon Disulfide, saturated atmosphere  2 days 17.4 11.8 Warped  *Chlorine, dry 1 week 8.8 2.2 Crazed and brittle  *Chlorine, moist 1 week 7.8 0.1 Crazed and brittle  *Chlorine, saturated solution 1 week — — Considerably softened and swollen  *I,4-Dioxane — — Dissolved  *Ethylene Oxide, gas 10 minutes (41°C, 105°F) — Unchanged  *Ethylene Oxide, gas 1 day 20.9 25.6 Swollen and softened  *Formaldehyde, 4% 10 min per day 0.2 — Unchanged for 5 days  *Formaldehyde, 35% 2 months 13.0 6.7 Swollen and softened  *Furfural — — Dissolved  Hydrogen Peroxide, 3% 1 year 1.7 1.1 Unchanged  Hydrogen Peroxide, 5% 2 days 1.4 1.3 Unchanged  Hydrogen Sulfide, dry 2 months 3.0 1.7 Unchanged  Hydrogen Sulfide, saturated solution 2 months 5.9 2.3 Unchanged  Hydroquinone, 20 g per gallon 1 week 2.4 1.0 Slightly stained yellow  *Methyl Methacrylate Monomer — — Dissolved  Ozone, 0.67-0.15 ppm 45 days (outdoors) — Unchanged  *Sulfur Dioxide, dry 2 months 19.4 8.6 Swollen, slightly warped  *Sulfur Dioxide, moist 2 months 31.9 10.2 Considerably swollen and warped  *Sulfur Dioxide, saturated 2 months 23.2 18.1	_		_	_		
•Carbon Disulfide, 1 week 25.8 1.6 •Carbon Disulfide, saturated atmosphere  2 days 17.4 11.8  •Chlorine, dry 1 week 8.8 2.2 Crazed and brittle •Chlorine, moist 1 week 7.8 0.1 Crazed and brittle •Chlorine, saturated solution 1 week — Considerably softened and swollen  •1,4-Dioxane — Dissolved  Ethylene Oxide, gas  10 minutes (41°C, 105°F) — Unchanged •Ethylene Oxide, gas 1 day 20.9 25.6  Formaldehyde, 4% 10 min per day 0.2 — Unchanged for 5 days  •Formaldehyde, 35% 2 months 13.0 6.7  •Furfural — Dissolved  Hydrogen Peroxide, 3% 1 year 1.7 1.1 Unchanged  Hydrogen Sulfide, dry 2 months 2.3 0.8 Unchanged  Hydrogen Sulfide, moist 2 months 3.0 1.7 Unchanged  Hydroquinone, 20 g per gallon 1 week 2.4 1.0  •Methyl Methacrylate Monomer — Dissolved  Ozone, 0.05-0.15 ppm 45 days (outdoors)— Ozone, 0.07-0pm 45 days (49°C, 120°F)— Yellowed  •Phenol 1 week — Dissolved  Sulfur, solid 1 week — Dissolved  •Sulfur Dioxide, moist 2 months 31.9 10.2  •Sulfur Dioxide, saturated 2 months 23.2 18.1		aid 6 months	<b>19.3</b>			
•Carbon Disulfide, saturated atmosphere			2.7			
• Chlorine, dry         1 week         8.8         2.2         Crazed and brittle           • Chlorine, moist         1 week         7.8         0.1         Crazed and brittle           • Chlorine, saturated solution 1 week         —         Considerably softened and swollen           • Italian         —         Dissolved           • Italian         —         Dissolved           • Ethylene Oxide, gas         —         —         Unchanged           • Ethylene Oxide, gas 1 day         20.9         25.6         Swollen and softened           Formaldehyde, 35% 2 months         13.0         6.7         Swollen and softened           • Formaldehyde, 35% 2 months         13.0         6.7         Swollen and softened           • Formaldehyde, 35% 2 months         13.0         6.7         Swollen and softened           • Formaldehyde, 35% 2 days         1.4         1.3         Unchanged           • Hydrogen Peroxide, 3% 1 year         1.7         1.1         Unchanged           Hydrogen Sulfide, dry 2 months         2.3         0.8         Unchanged           Hydrogen Sulfide, saturated solution         2         2.3         Unchanged           Hydroquinone, 20 g per gallon 1 week 2.4 1.0         Slightly stained yellow         • Sulfurly Methacrylate Monomer<		2 11 0 0 2 2			Softened and swollen	
•Chlorine, dry 1 week 8.8 2.2 Crazed and brittle •Chlorine, moist 1 week 7.8 0.1 Crazed and brittle •Chlorine, saturated solution 1 week — — Considerably softened and swollen  •1,4-Dioxane — — Dissolved  •thylene Oxide, gas  10 minutes (41°C, 105°F) — — Unchanged  •Ethylene Oxide, gas 1 day 20.9 25.6 Swollen and softened  Formaldehyde, 4% 10 min per day 0.2 — Unchanged for 5 days  •Formaldehyde, 35% 2 months 13.0 6.7 Swollen and softened  •Furfural — — Dissolved  Hydrogen Peroxide, 3% 1 year 1.7 1.1 Unchanged  Hydrogen Sulfide, dry 2 months 2.3 0.8 Unchanged  Hydrogen Sulfide, moist 2 months 3.0 1.7 Unchanged  Hydrogen Sulfide, saturated solution  2 months 5.9 2.3 Unchanged  Hydroquinone, 20 g per gallon 1 week 2.4 1.0 Slightly stained yellow  •Methyl Methacrylate Monomer — — Dissolved  Ozone, 0.05-0.15 ppm 45 days (outdoors)—— Unchanged  •Sulfur, solid 1 week — — Decomposed  •Sulfur, solid 1 week — — Decomposed  •Sulfur Dioxide, moist 2 months 31.9 10.2  •Sulfur Dioxide, saturated 2 months 23.2 18.1  •Swollen and brittle  Carzed and brittle  Carzed and brittle  Carzed and brittle  Carzed and brittle  Considerably softened and swollen  Bissolved  Unchanged  Value 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1 1	<ul> <li>Carbon Disulfide,</li> </ul>		ospher	$\mathbf{e}$		
•Chlorine, moist 1 week 7.8 0.1 Crazed and brittle •Chlorine, saturated solution 1 week — — Considerably softened and swollen •1,4-Dioxane — — Dissolved  Ethylene Oxide, gas  10 minutes (41°C, 105°F) — — Unchanged •Ethylene Oxide, gas 1 day 20.9 25.6 Swollen and softened Formaldehyde, 4% 10 min per day 0.2 — Unchanged for 5 days •Formaldehyde, 35% 2 months 13.0 6.7 Swollen and softened  Formaldehyde, 35% 2 months 1.7 1.1 Unchanged  Hydrogen Peroxide, 5% 2 days 1.4 1.3 Unchanged Hydrogen Peroxide, 5% 2 days 1.4 1.3 Unchanged Hydrogen Sulfide, dry 2 months 2.3 0.8 Unchanged Hydrogen Sulfide, saturated solution  2 months 5.9 2.3 Unchanged Hydroquinone, 20 g per gallon 1 week 2.4 1.0 Slightly stained yellow •Methyl Methacrylate Monomer — — Dissolved  Ozone, 0.05-0.15 ppm 45 days (outdoors) — Unchanged  •Nitrobenzene — — Dissolved  Ozone, 0.7 ppm 45 days (49°C, 120°F) — Yellowed •Phenol 1 week — — Decomposed •Stlyrene Monomer — — Dissolved  Sulfur, solid 1 week — — Decomposed •Sulfur Dioxide, dry 2 months 19.4 8.6 Swollen, slightly warped  •Sulfur Dioxide, moist 2 months 31.9 10.2  •Sulfur Dioxide, saturated 2 months 23.2 18.1  Crazed and brittle Considerably softened and swollen  Sulfur Dioxide, saturated 2 months 23.2 18.1		$2  ext{ days}$	<b>17.4</b>	11.8		
•Chlorine, saturated solution 1 week — — Considerably softened and swollen  •1,4-Dioxane — — Dissolved  Ethylene Oxide, gas  10 minutes (41°C, 105°F) — — Unchanged  •Ethylene Oxide, gas 1 day 20.9 25.6 Swollen and softened  Formaldehyde, 4% 10 min per day 0.2 — Unchanged for 5 days  •Formaldehyde, 35% 2 months 13.0 6.7 Swollen and softened  •Furfural — — Dissolved  Hydrogen Peroxide, 5% 2 days 1.4 1.3 Unchanged  Hydrogen Peroxide, 5% 2 days 1.4 1.3 Unchanged  Hydrogen Sulfide, dry 2 months 2.3 0.8 Unchanged  Hydrogen Sulfide, saturated solution  2 months 5.9 2.3 Unchanged  Hydroquinone, 20 g per gallon 1 week 2.4 1.0 Slightly stained yellow  •Methyl Methacrylate Monomer — — Dissolved  Ozone, 0.05-0.15 ppm 45 days (outdoors)—— Unchanged  •Nitrobenzene — — Dissolved  Ozone, 0.7 ppm 45 days (49°C, 120°F)— — Yellowed  •Phenol 1 week — — Decomposed  •Styrene Monomer — — Dissolved  •Sulfur, solid 1 week — — Decomposed  •Sulfur Dioxide, dry 2 months 31.9 10.2 Considerably swollen and warped  •Sulfur Dioxide, saturated 2 months 23.2 18.1 Swollen and warped		1  week	8.8	2.2	Crazed and brittle	
**Swollen**  **I,4-Dioxane**  **I0 minutes (41°C, 105°F) — Unchanged  **Ethylene Oxide, gas 1 day 20.9 25.6  **Ethylene Oxide, gas 1 day 20.9 25.6  **Formaldehyde, 4% 10 min per day 0.2 — Unchanged for 5 days  **Formaldehyde, 35% 2 months 13.0 6.7  **Furfural — — Dissolved  Hydrogen Peroxide, 3% 1 year 1.7 1.1 Unchanged  Hydrogen Peroxide, 5% 2 days 1.4 1.3 Unchanged  Hydrogen Sulfide, dry 2 months 2.3 0.8 Unchanged  Hydrogen Sulfide, moist 2 months 3.0 1.7 Unchanged  Hydroquinone, 20 g per gallon 1 week 2.4 1.0 Slightly stained yellow  **Methyl Methacrylate Monomer — — Dissolved  Ozone, 0.05-0.15 ppm 45 days (outdoors)—— Unchanged  **Purfured — Dissolved  Ozone, 0.7 ppm 45 days (49°C, 120°F)— — Yellowed  **Styrene Monomer — — Dissolved  **Sulfur, solid 1 week — — Decomposed  **Sulfur, solid 1 week — — Unchanged  **Sulfur Dioxide, dry 2 months 19.4 8.6 Swollen, slightly warped  **Sulfur Dioxide, moist 2 months 31.9 10.2 Considerably swollen and warped  **Sulfur Dioxide, saturated 2 months 23.2 18.1 Swollen and warped	The state of the s			0.1		
**Hydrogen Peroxide, gas 1 day 20.9 25.6 Swollen and softened Formaldehyde, 4% 10 min per day 0.2 — Unchanged of Formaldehyde, 35% 2 months 13.0 6.7 Swollen and softened Formaldehyde, 35% 2 months 13.0 6.7 Swollen and softened Hydrogen Peroxide, 3% 1 year 1.7 1.1 Unchanged Hydrogen Peroxide, 5% 2 days 1.4 1.3 Unchanged Hydrogen Sulfide, dry 2 months 2.3 0.8 Unchanged Hydrogen Sulfide, moist 2 months 3.0 1.7 Unchanged Hydrogen Sulfide, saturated solution 2 months 5.9 2.3 Unchanged Hydroquinone, 20 g per gallon 1 week 2.4 1.0 Slightly stained yellow •Methyl Methacrylate Monomer — Dissolved Ozone, 0.05-0.15 ppm 45 days (outdoors)— Unchanged •Nitrobenzene — Dissolved •Phenol 1 week — Decomposed •Styrene Monomer — Dissolved •Styrene Monomer — Dissolved •Sulfur, solid 1 week — Unchanged •Sulfur Dioxide, dry 2 months 19.4 8.6 Swollen, slightly warped •Sulfur Dioxide, moist 2 months 23.2 18.1 Swollen and warped	•Chlorine, saturate	ed solution 1 we	eek —	_	Considerably softened and	
Ethylene Oxide, gas  10 minutes (41°C, 105°F) — Unchanged  •Ethylene Oxide, gas 1 day  20.9  •Ethylene Oxide, gas 1 day  20.9  •Ethylene Oxide, gas 1 day  •Ethylene Oxide, gas 1 day  20.9  25.6  Swollen and softened  Unchanged for 5 days  •Formaldehyde, 35% 2 months  •Furfural  ———————————————————————————————————						
*Ethylene Oxide, gas 1 day 20.9 25.6 Swollen and softened Formaldehyde, 4% 10 min per day 0.2 — Unchanged for 5 days Formaldehyde, 35% 2 months 13.0 6.7 Swollen and softened Furfural — Dissolved Hydrogen Peroxide, 5% 2 days 1.4 1.3 Unchanged Hydrogen Peroxide, 5% 2 days 1.4 1.3 Unchanged Hydrogen Sulfide, dry 2 months 2.3 0.8 Unchanged Hydrogen Sulfide, moist 2 months 3.0 1.7 Unchanged Hydrogen Sulfide, saturated solution 2 months 5.9 2.3 Unchanged Hydroquinone, 20 g per gallon 1 week 2.4 1.0 Slightly stained yellow Methyl Methacrylate Monomer — Dissolved Ozone, 0.05-0.15 ppm 45 days (outdoors)—— Unchanged Ozone, 0.7 ppm 45 days (49°C, 120°F)— Yellowed Phenol 1 week — Decomposed Sulfur, solid 1 week — Dissolved Unchanged Sulfur, solid 1 week — Unchanged Sulfur Dioxide, dry 2 months 19.4 8.6 Swollen, slightly warped Sulfur Dioxide, moist 2 months 23.2 18.1 Swollen and warped	· ·		_	_	Dissolved	
•Ethylene Oxide, gas 1 day Formaldehyde, 4% 10 min per day 0.2  •Formaldehyde, 35% 2 months  •Formaldehyde, 35% 2 months  •Furfural						
Formaldehyde, 4% 10 min per day 0.2 — Unchanged for 5 days  Formaldehyde, 35% 2 months 13.0 6.7 Swollen and softened  Furfural — — Dissolved  Hydrogen Peroxide, 3% 1 year 1.7 1.1 Unchanged  Hydrogen Peroxide, 5% 2 days 1.4 1.3 Unchanged  Hydrogen Sulfide, dry 2 months 2.3 0.8 Unchanged  Hydrogen Sulfide, moist 2 months 3.0 1.7 Unchanged  Hydrogen Sulfide, saturated solution  2 months 5.9 2.3 Unchanged  Hydroquinone, 20 g per gallon 1 week 2.4 1.0 Slightly stained yellow  Methyl Methacrylate Monomer — — Dissolved  Ozone, 0.05-0.15 ppm 45 days (outdoors)—— Unchanged  Ozone, 0.7 ppm 45 days (49°C, 120°F)— — Yellowed  Phenol 1 week — — Decomposed  Styrene Monomer — — Dissolved  Sulfur, solid 1 week — — Unchanged  Sulfur Dioxide, dry 2 months 19.4 8.6 Swollen, slightly warped  Sulfur Dioxide, moist 2 months 31.9 10.2 Considerably swollen and warped  Sulfur Dioxide, saturated 2 months 23.2 18.1					_	
•Formaldehyde, 35% 2 months  •Furfural  ———————————————————————————————————		•		<b>25.6</b>		
•Furfural — — — Dissolved Hydrogen Peroxide, 3% 1 year 1.7 1.1 Unchanged Hydrogen Peroxide, 5% 2 days 1.4 1.3 Unchanged Hydrogen Sulfide, dry 2 months 2.3 0.8 Unchanged Hydrogen Sulfide, moist 2 months 3.0 1.7 Unchanged Hydrogen Sulfide, saturated solution  2 months 5.9 2.3 Unchanged Hydroquinone, 20 g per gallon 1 week 2.4 1.0 Slightly stained yellow •Methyl Methacrylate Monomer — — Dissolved •Nitrobenzene — — Dissolved Ozone, 0.05-0.15 ppm 45 days (outdoors)—— Unchanged Ozone, 0.7 ppm 45 days (49°C, 120°F)— — Yellowed •Phenol 1 week — — Decomposed •Styrene Monomer — — Dissolved Sulfur, solid 1 week — — Dissolved •Sulfur Dioxide, dry 2 months 19.4 8.6 Swollen, slightly warped •Sulfur Dioxide, moist 2 months 31.9 10.2 Considerably swollen and warped •Sulfur Dioxide, saturated 2 months 23.2 18.1						
Hydrogen Peroxide, 3% 1 year 1.7 1.1 Unchanged Hydrogen Peroxide, 5% 2 days 1.4 1.3 Unchanged Hydrogen Sulfide, dry 2 months 2.3 0.8 Unchanged Hydrogen Sulfide, moist 2 months 3.0 1.7 Unchanged Hydrogen Sulfide, saturated solution  2 months 5.9 2.3 Unchanged Hydroquinone, 20 g per gallon 1 week 2.4 1.0 Slightly stained yellow •Methyl Methacrylate Monomer — — Dissolved •Nitrobenzene — — Dissolved Ozone, 0.05-0.15 ppm 45 days (outdoors)—— Unchanged Ozone, 0.7 ppm 45 days (49°C, 120°F)— — Yellowed •Phenol 1 week — — Decomposed •Styrene Monomer — — Dissolved Sulfur, solid 1 week — — Dissolved Sulfur Dioxide, dry 2 months 19.4 8.6 Swollen, slightly warped •Sulfur Dioxide, moist 2 months 31.9 10.2 Considerably swollen and warped •Sulfur Dioxide, saturated 2 months 23.2 18.1		% 2 months	13.0	<b>6.7</b>		
Hydrogen Peroxide, 5% 2 days 1.4 1.3 Unchanged Hydrogen Sulfide, dry 2 months 2.3 0.8 Unchanged Hydrogen Sulfide, moist 2 months 3.0 1.7 Unchanged Hydrogen Sulfide, saturated solution  2 months 5.9 2.3 Unchanged Hydroquinone, 20 g per gallon 1 week 2.4 1.0 Slightly stained yellow  • Methyl Methacrylate Monomer — — Dissolved • Nitrobenzene — — Dissolved Ozone, 0.05-0.15 ppm 45 days (outdoors)—— Unchanged Ozone, 0.7 ppm 45 days (49°C, 120°F)— — Yellowed • Phenol 1 week — — Decomposed • Styrene Monomer — — Dissolved Sulfur, solid 1 week — — Unchanged • Sulfur Dioxide, dry 2 months 19.4 8.6 Swollen, slightly warped • Sulfur Dioxide, moist 2 months 31.9 10.2 Considerably swollen and warped • Sulfur Dioxide, saturated 2 months 23.2 18.1 Swollen and warped			_	_		
Hydrogen Sulfide, dry 2 months 2.3 0.8 Unchanged Hydrogen Sulfide, moist 2 months 3.0 1.7 Unchanged Hydrogen Sulfide, saturated solution  2 months 5.9 2.3 Unchanged Hydroquinone, 20 g per gallon 1 week 2.4 1.0 Slightly stained yellow  • Methyl Methacrylate Monomer — — Dissolved • Nitrobenzene — — Dissolved Ozone, 0.05-0.15 ppm 45 days (outdoors)—— Unchanged Ozone, 0.7 ppm 45 days (49°C, 120°F)— — Yellowed • Phenol 1 week — — Decomposed • Styrene Monomer — — Dissolved Sulfur, solid 1 week — — Unchanged • Sulfur Dioxide, dry 2 months 19.4 8.6 Swollen, slightly warped • Sulfur Dioxide, moist 2 months 31.9 10.2 Considerably swollen and warped • Sulfur Dioxide, saturated 2 months 23.2 18.1 Swollen and warped	•				_	
Hydrogen Sulfide, moist 2 months 3.0 1.7 Unchanged  Hydrogen Sulfide, saturated solution  2 months 5.9 2.3 Unchanged  Hydroquinone, 20 g per gallon 1 week 2.4 1.0 Slightly stained yellow  •Methyl Methacrylate Monomer — — Dissolved  •Nitrobenzene — — Dissolved  Ozone, 0.05-0.15 ppm 45 days (outdoors)—— Unchanged  Ozone, 0.7 ppm 45 days (49°C, 120°F)— — Yellowed  •Phenol 1 week — — Decomposed  •Styrene Monomer — Dissolved  Sulfur, solid 1 week — — Unchanged  •Sulfur Dioxide, dry 2 months 19.4 8.6 Swollen, slightly warped  •Sulfur Dioxide, moist 2 months 31.9 10.2 Considerably swollen and warped  •Sulfur Dioxide, saturated 2 months 23.2 18.1 Swollen and warped	-			<u> </u>		
Hydrogen Sulfide, saturated solution  2 months 5.9 2.3 Unchanged  Hydroquinone, 20 g per gallon 1 week 2.4 1.0 Slightly stained yellow  •Methyl Methacrylate Monomer — — Dissolved  •Nitrobenzene — — Dissolved  Ozone, 0.05-0.15 ppm 45 days (outdoors)—— Unchanged  Ozone, 0.7 ppm 45 days (49°C, 120°F)— — Yellowed  •Phenol 1 week — — Decomposed  •Styrene Monomer — Dissolved  Sulfur, solid 1 week — — Unchanged  •Sulfur Dioxide, dry 2 months 19.4 8.6 Swollen, slightly warped  •Sulfur Dioxide, moist 2 months 31.9 10.2 Considerably swollen and warped  •Sulfur Dioxide, saturated 2 months 23.2 18.1 Swollen and warped					<u> </u>	
2 months 5.9 2.3 Unchanged Hydroquinone, 20 g per gallon 1 week 2.4 1.0 Slightly stained yellow  •Methyl Methacrylate Monomer — — Dissolved  •Nitrobenzene — — Dissolved Ozone, 0.05-0.15 ppm 45 days (outdoors)—— Unchanged Ozone, 0.7 ppm 45 days (49°C, 120°F)— — Yellowed  •Phenol 1 week — — Decomposed  •Styrene Monomer — — Dissolved Sulfur, solid 1 week — — Unchanged  •Sulfur Dioxide, dry 2 months 19.4 8.6 Swollen, slightly warped  •Sulfur Dioxide, moist 2 months 31.9 10.2 Considerably swollen and warped  •Sulfur Dioxide, saturated 2 months 23.2 18.1 Swollen and warped	• •			1.7	Unchanged	
Hydroquinone, 20 g per gallon 1 week 2.4 1.0  •Methyl Methacrylate Monomer — — Dissolved  •Nitrobenzene — — Dissolved Ozone, 0.05-0.15 ppm 45 days (outdoors)—— Unchanged Ozone, 0.7 ppm 45 days (49°C, 120°F)— — Yellowed  •Phenol 1 week — — Decomposed  •Styrene Monomer — — Dissolved Sulfur, solid 1 week — — Unchanged  •Sulfur Dioxide, dry 2 months 19.4 8.6 Swollen, slightly warped  •Sulfur Dioxide, moist 2 months 31.9 10.2 Considerably swollen and warped  •Sulfur Dioxide, saturated 2 months 23.2 18.1 Swollen and warped						
<ul> <li>•Methyl Methacrylate Monomer — — Dissolved</li> <li>•Nitrobenzene — — Dissolved</li> <li>Ozone, 0.05-0.15 ppm 45 days (outdoors)—— Unchanged</li> <li>Ozone, 0.7 ppm 45 days (49°C, 120°F)— — Yellowed</li> <li>•Phenol 1 week — — Decomposed</li> <li>•Styrene Monomer — — Dissolved</li> <li>Sulfur, solid 1 week — — Unchanged</li> <li>•Sulfur Dioxide, dry 2 months 19.4 8.6 Swollen, slightly warped</li> <li>•Sulfur Dioxide, moist 2 months 31.9 10.2 Considerably swollen and warped</li> <li>•Sulfur Dioxide, saturated 2 months 23.2 18.1</li> <li>Swollen and warped</li> </ul>					_	
<ul> <li>Nitrobenzene — — Dissolved</li> <li>Ozone, 0.05-0.15 ppm 45 days (outdoors)—— Unchanged</li> <li>Ozone, 0.7 ppm 45 days (49°C, 120°F)— — Yellowed</li> <li>Phenol 1 week — — Decomposed</li> <li>Styrene Monomer — Dissolved</li> <li>Sulfur, solid 1 week — — Unchanged</li> <li>Sulfur Dioxide, dry 2 months 19.4 8.6 Swollen, slightly warped</li> <li>Sulfur Dioxide, moist 2 months 31.9 10.2 Considerably swollen and warped</li> <li>Sulfur Dioxide, saturated 2 months 23.2 18.1</li> <li>Swollen and warped</li> </ul>			veek 2.4	4 1.0		
Ozone, 0.05-0.15 ppm 45 days (outdoors)—— Unchanged Ozone, 0.7 ppm 45 days (49°C, 120°F)— — Yellowed  •Phenol 1 week — — Decomposed •Styrene Monomer — — Dissolved Sulfur, solid 1 week — — Unchanged •Sulfur Dioxide, dry 2 months 19.4 8.6 Swollen, slightly warped •Sulfur Dioxide, moist 2 months 31.9 10.2 Considerably swollen and warped •Sulfur Dioxide, saturated 2 months 23.2 18.1 Swollen and warped	-	ate Monomer	_	_		
Ozone, 0.7 ppm 45 days (49°C, 120°F)— — Yellowed  •Phenol 1 week — — Decomposed  •Styrene Monomer — — Dissolved Sulfur, solid 1 week — — Unchanged  •Sulfur Dioxide, dry 2 months 19.4 8.6 Swollen, slightly warped  •Sulfur Dioxide, moist 2 months 31.9 10.2 Considerably swollen and warped  •Sulfur Dioxide, saturated 2 months 23.2 18.1 Swollen and warped			_ \	_		
<ul> <li>Phenol 1 week — — Decomposed</li> <li>Styrene Monomer — Dissolved</li> <li>Sulfur, solid 1 week — — Unchanged</li> <li>Sulfur Dioxide, dry 2 months 19.4 8.6 Swollen, slightly warped</li> <li>Sulfur Dioxide, moist 2 months 31.9 10.2 Considerably swollen and warped</li> <li>Sulfur Dioxide, saturated 2 months 23.2 18.1 Swollen and warped</li> </ul>	,	•			_	
•Styrene Monomer — — Dissolved Sulfur, solid 1 week — — Unchanged •Sulfur Dioxide, dry 2 months 19.4 8.6 Swollen, slightly warped •Sulfur Dioxide, moist 2 months 31.9 10.2 Considerably swollen and warped •Sulfur Dioxide, saturated 2 months 23.2 18.1 Swollen and warped		•	'F)—	_		
Sulfur, solid 1 week — — Unchanged  •Sulfur Dioxide, dry 2 months 19.4 8.6 Swollen, slightly warped  •Sulfur Dioxide, moist 2 months 31.9 10.2 Considerably swollen and warped  •Sulfur Dioxide, saturated 2 months 23.2 18.1 Swollen and warped			_	_		
•Sulfur Dioxide, dry 2 months 19.4 8.6 Swollen, slightly warped Considerably swollen and warped •Sulfur Dioxide, saturated 2 months 23.2 18.1 Swollen and warped			_	_		
•Sulfur Dioxide, moist 2 months 31.9 10.2 Considerably swollen and warped •Sulfur Dioxide, saturated 2 months 23.2 18.1 Swollen and warped			<del>_</del> .			
•Sulfur Dioxide, saturated 2 months 23.2 18.1 Swollen and warped		v 2 months	<b>19.4</b>	<b>8.6</b>	Swollen, slightly warped	
•Sulfur Dioxide, saturated 2 months 23.2 18.1 Swollen and warped						
	<del>-</del>					
Chemical Compatibility Chart continued		oist 2 months			warped	
		oist 2 months turated 2 mont	ths 23.2	2 18.1	warped	

**Percent Increase Observed Condition** Time Exposed Weight Thickness of Plastic •Sulfur Dioxide in Hydrocarbons — Dissolved •Sulfur Dioxide and Hydrocarbon Vapor 2 months 19.2 Swollen 11.5 •Titanium Tetrachloride 3 days Very brittle Trinitrotoluene (TNT), water slurry 4 weeks Stained

### **Optional Accessories**

Your phantom may include some of the following accessories:

#### **Tumor Rotation Ball**

Some of the Liqui-Phil<sup>TM</sup> Phantoms that we manufacture have rotation ball tumor supports. For proper functioning, it is important that the internal o-ring within the rotation ball is lightly lubricated with petroleum jelly or silicon lubricant. The tumor and organ support rods will fit through the rotation balls and allow for varied positioning. The rotation ball lock nut holds the tumor into place and also seals the phantom in order to prevent leaking.

To insert tumors into the organs or phantom, remove the port plug and assemble the tumor through the rotation ball port and tumor into the organ. Not all tumor sizes will fit into all organs. When organs or tumors are not required in a study, the port can be sealed by rotating the rotation ball so the internal hole is perpendicular to the port hole.

#### **Optional Tumor Vessels**

The tumor vessels come in many different shapes, sizes, and materials. It is recommended that a test be performed prior to filling with solutions. This should be done by applying a small amount of the desired solution onto the outside of the tumor vessel and observing the material for any changes in appearance. If unsure of the solution's affect on the plastic, it is recommended that the tumor vessels be emptied and cleaned after each use.