

part number
SK150/100

page 1 of 1



Sectional Head Phantoms, SK150 and SK100

The Phantom Laboratory manufactures two Sectional Head Phantoms. The SK150 is manufactured with a natural human skull and upper cervical vertebrae. The SK100 is manufactured with a natural human skull.

The natural bones are cast from our RANDO® material with an internal air cavity representing the oral, pharynx and trachea anatomy. The opaque RANDO® material simulates the x-ray absorbency, atomic number and specific gravity of soft human tissue. An additional benefit of this urethane material is its durability. The Sectional Phantoms will withstand substantial impacts and continuous handling without damage.



The Phantom Laboratory

phone
800-525-1190 or 518-692-1190
fax
518-692-3329

mail
P.O. Box 511
Salem, New York 12865
shipping
2727 State Route 29
Greenwich, New York 12834

email
info@phantomlab.com
web
www.phantomlab.com