

part number  
XA231

page 1 of 1



## Sectional Hand Phantoms, XA231R and XA231L

There are two Sectional Hand Phantoms available: XA231L, left hand pronated, and XA231R, right hand in relaxed position.

For each Sectional Phantom a set of natural human bones are reconstructed into an average size male mold and then cast with our RANDO® material. The opaque RANDO® material simulates the x-ray absorbency, atomic number and specific gravity of soft human tissue. An additional benefit of this urethane material is its durability. The Sectional Phantoms will withstand substantial impacts and continuous handling without damage.

Many human skeletons are not the same size and shape as our molds and reflect natural human characteristics such as lack of symmetry and distorted joints. As our technicians reconstruct the skeleton, minor adjustments may be made to facilitate positioning within the mold.



## The Phantom Laboratory

**phone**  
800-525-1190 or 518-692-1190  
**fax**  
518-692-3329

**mail**  
P.O. Box 511  
Salem, New York 12865  
**shipping**  
2727 State Route 29  
Greenwich, New York 12834

**email**  
info@phantomlab.com  
**web**  
www.phantomlab.com