


part number
CTP299

page 1 of 1



Teflon Annulus, CTP299

The Teflon annulus can be slid over a 20cm Catphan® housing to duplicate the beam hardening effects created by the skull in head imaging. The Teflon annulus is approximately 7mm thick.



The Phantom Laboratory

phone
800-525-1190 or 518-692-1190
fax
518-692-3329

mail
P.O. Box 511
Salem, New York 12865
shipping
2727 State Route 29
Greenwich, New York 12834

email
info@phantomlab.com
web
www.phantomlab.com