


**part number**  
**CTP755**


page 1 of 1



### Scatter Extension for Catphan® 600

The CTP755 scatter extension adds 5 cm of extra uniformity material on the end of the Catphan® 600. This material helps account for the scatter effects seen with wide multi slice beam or cone beam CT protocols. The scatter extension will have minimal or no effect when the rotated beam covers less than 40 mm in the z direction.

The CTP755 only mounts on the Catphan® 600 and cannot be used with other Catphan® models.



### The Phantom Laboratory

**phone**  
800-525-1190 or 518-692-1190  
**fax**  
518-692-3329

**mail**  
P.O. Box 511  
Salem, New York 12865  
**shipping**  
2727 State Route 29  
Greenwich, New York 12834

**email**  
info@phantomlab.com  
**web**  
www.phantomlab.com