


part number
See Table

page 1 of 1



Catphan® Body Annuli

Catphan® annuli are cast from uniformity material. Annuli are available in a variety of round and oval sizes and mimic the attenuation of patients and their effect on image quality. All annuli are designed to be slipped over 20cm OD Catphan® phantoms.

Part Number	Shape	H (cm)	W (cm)	L (cm)
CTP539	round	30	30	5
CTP653	round	32	32	5
CTP653-10	round	32	32	10
CTP653-15	round	32	32	15
CTP653-20	round	32	32	20
CTP540	round	35	35	5
CTP650	round	38	38	5
CTP650-10	round	38	38	10
CTP650-15	round	38	38	15
CTP650-20	round	38	38	20
CTP579	oval	25	35	5
CTP579-10	oval	25	35	10
CTP579-15	oval	25	35	15
CTP579-20	oval	25	35	20
CTP651	oval	30	38	5
CTP651-10	oval	30	38	10
CTP651-15	oval	30	38	15
CTP651-20	oval	30	38	20
CTP750-10	oval	42	55	10
CTP750-20	oval	42	55	20



The Phantom Laboratory

phone
800-525-1190 or 518-692-1190
fax
518-692-3329

mail
P.O. Box 511
Salem, New York 12865
shipping
2727 State Route 29
Greenwich, New York 12834

email
info@phantomlab.com
web
www.phantomlab.com